

88000  
to  
88040

## SAFESEAL® QUATTRO STERILIZATION POUCHES Class I, Medical Device

### FEATURES & BENEFITS

- For use as packing material for medical devices in sterilization by steam, ethylene oxide gas in healthcare establishments
- Non-tearing, thick multi-layer film made of a combination of polypropylene and polyethylene.
- Transparent, blue tinted film facilitates both easy instrument viewing and identification of puncture or tears in the film
- Clean, fibre-free separation of paper and film ensures safe presentation of sterile products.
- Pouches have a wide adhesive sealing tape which easily and securely seals the pouch providing a reliable seal until pouch is opened.
- Permit to keep the sterilization conditions of the inner product after sterilization : up to 1 year.



Process Indicator Colour Change Guide :



Steam: changing colour from pink to cocoa



EO: changing colour from blue to golden brown

### PRODUCT DESCRIPTION

Product name: **SafeSeal® Quattro Sterilization Pouches**  
 Product type: Single use  
 Material: Paper, polypropylene and polyethylene  
 Medical grade paper: 60 g/m<sup>2</sup>  
 Film : 50µm

Sterilization Cycle Parameter	Gravity Steam		Ethylene Oxide Gas
	30	4	
Time (min)	30	4	100-120
Temperature (°C)	121	134	50
Relative humidity (%)	N/A	N/A	60-85
Sterilant Concentration	N/A	N/A	600
Drying-Aeration time	15-30 min	20-30 min	12 hrs

N/A = Not applicable

For color change guide, facts on the Stated Value (SV) and all other information, consult [www.medicom.com](http://www.medicom.com)

Item	Size mm
88000	70 x 229
88005	89 x 133
88010	89 x 229
88015	57 x 102
88025	133 x 254
88030	190 x 330
88035	254 x 356
88040	305 x 432



### INDUSTRIES

Medical / Dental

88000  
to  
88040

# SAFESEAL® QUATTRO STERILIZATION POUCHES

Class I, Medical Device

## CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC  
Manufactured under ISO 13485 certified quality system

In conformity with the relevant harmonized standards EN 868-5, EN ISO 11607-1 and international Standard ISO 11140-1

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 36 months

## LOGISTIC INFORMATION



### Specifications outer case

Item	Size mm	Gross weight Kg	QTY/ pallet
88000	268 x 160 x 202	2,9	180
88005	483 x 173 x 74	2,1	200
88010	298 x 267 x 138	3,7	120
88015	328 x 141 x 138	2,2	285
88025	298 x 298 x 352	10,0	48
88030	328 x 210 x 370	8,5	60
88035	394 x 330 x 267	12,3	42
88040	476 x 337 x 343	17,7	25



### Specifications inner case

Size mm	QTY
260 x 64 x 76	6 x 200
165 x 64 x 95	5 x 200
260 x 64 x 95	6 x 200
133 x 64 x 64	10 x 200
286 x 64 x 140	10 x 200
362 x 64 x 197	5 x 200
387 x 64 x 260	5 x 200
464 x 64 x 311	5 x 200



### EAN Codes

Item	Inner case	Outer case
88000	8 718 007 900 422	8 718 007 900 415
88005	8 718 007 900 446	8 718 007 900 439
88010	8 718 007 900 569	8 718 007 900 552
88015	8 718 007 900 811	8 718 007 900 576
88025	8 718 007 900 835	8 718 007 900 828
88030	8 718 007 900 859	8 718 007 900 842
88035	8 718 007 900 873	8 718 007 900 866
88040	8 718 007 900 903	8 718 007 900 880



KOLMI HOPEN SASU

Bld de la chanterie | 49124 Saint Barthélemy d'Anjou – France

Bld de la chanterie | BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex – France

T : + 33 (0) 241 963 434 | F : + 33 (0) 241 963 453 | sales@kolmi-hopen.com | www.h-open.com