


PROTAPER • GOLD™



The same winning technique with **greater flexibility**



THEN: PROTAPER® UNIVERSAL sets a new standard in efficiency



NOW: The PROTAPER® legacy continues with PROTAPER GOLD™

SAME: File Sequence • File Sizes • Motor Settings

DIFFERENT: Heat Treatment • Gold Colour

IMPROVED: Flexibility • Cyclic Fatigue Resistance • Tactile Feeling • Safety

**+
WE
KNOW
ENDO.**

**PROTAPER GOLD™
IS NOW REPLACING
PROTAPER® UNIVERSAL**

PROTAPER•GOLD™

Product references

PROTAPER GOLD™ FILES



SX STERILE
19MM A 0410 219 G01 03



S1 STERILE
21MM A 0410 221 G01 03
25MM A 0410 225 G01 03
31MM A 0410 231 G01 03



F3 STERILE
21MM A 0411 221 G03 03
25MM A 0411 225 G03 03
31MM A 0411 231 G03 03



S2 STERILE
21MM A 0410 221 G02 03
25MM A 0410 225 G02 03
31MM A 0410 231 G02 03



F4 STERILE
21MM A 0411 221 G04 03
25MM A 0411 225 G04 03
31MM A 0411 231 G04 03



F1 STERILE
21MM A 0411 221 G01 03
25MM A 0411 225 G01 03
31MM A 0411 231 G01 03



F5 STERILE
21MM A 0411 221 G05 03
25MM A 0411 225 G05 03
31MM A 0411 231 G05 03



F2 STERILE
21MM A 0411 221 G02 03
25MM A 0411 225 G02 03
31MM A 0411 231 G02 03

ASSORTMENT FILE PACKS
21MM SX/F3 A 0409 221 G01 03
25MM SX/F3 A 0409 225 G01 03
31MM SX/F3 A 0409 231 G01 03



“PROTAPER GOLD™ produces consistently predictable shapes for predictable cleaning and obturation that is easier and safer than ever before.”

Dr. John West²

All current PROTAPER® Universal obturation methods can be used with PROTAPER GOLD™.

PROTAPER GOLD™ now features an 11mm handle for improved intra oral visibility (versus 13mm with PROTAPER® Universal).

In some instances PROTAPER GOLD™ files may appear slightly curved. This feature is due to the GOLD technology allowing the files to be pre bent to achieve better access.



For more information, please visit www.dentsplymaillefer.com

1. Average value for PROTAPER GOLD™ compared to PROTAPER® Universal. Based on internal testing. Data on file.
2. PROTAPER GOLD™ (PTG) Shapers and Finishers. Dentsply Sirona (DTDS).