

**STRONGER, GLOSSIER AND MORE RELIABLE –  
A NEW GENERATION OF PROVEN QUALITY!**



**TEMPORARY CROWN AND BRIDGE MATERIAL**



## THE NEW GENERATION OF VISALYS® TEMP

**VISALYS® TEMP ONE:ONE** offers maximum stability, high fracture strength, and excellent edge stability. As a 1:1 crown and bridge material, it delivers precise, reliable, and aesthetically high-quality results – without compromise.



**Minimal waste, maximum efficiency. Ideal for modern practices with the highest standards.**

Only initial waste required  
**MULTIPLE APPLICATIONS  
POSSIBLE\***



Up to **62 ANTERIOR CROWNS**

**SAVE 35% MORE** on materials

**High  
fracture strength  
& edge stability.**

**Low  
smear layer:**  
Easier handling  
after removal.

**High gloss**  
without polishing.

**1:1  
SYSTEM**

**Precise milling:**  
Low dust generation,  
precise results.


**Easy to repair:**  
Easy to supplement  
after curing.




# RESTAURATION

... from **TEMPORARY** to **PERMANENT** – uncompromising strength

## OUTSTANDING FRACTURE STRENGTH



## CONVENIENT WORKING AND SETTING TIMES



**Note:** The working times apply to products stored and processed at a temperature of 23°C/74°F and 50 % relative humidity. Higher temperatures reduce the working time and lower temperatures extend these times.

**VISALYS® TEMP ONE:ONE**  
**GET TO KNOW**

# VISALYS® TEMP ONE:ONE – PRODUCT RANGE

Color	Normal pack	Bonus pack
	1 x 50 mL (73 g) cartridges, 15 orange mixing tips	5 x 50 mL (73 g) cartridges, 15 orange mixing tips
Bright (B1-A1)	REF 13850	
Medium (A2-A3)	REF 13851	REF 13856
Natural Bleach	REF 13849	
Bleach	REF 13853	
Mixing tips orange, 60 units		REF 17217



## EXTRAS



Silginat® Refill pack	2 x 380 mL cartridges	REF 14713
Dynamic mixers	45 units	REF 17900
Multi Tray Anterior	30 units	REF 17752
Multi Tray Posterior	50 units	REF 17750
Multi Tray Quadrant	30 units	REF 17753
Orthoskavident® C Normal pack	1 x 150 mL bottles	REF 13063
Orthoskavident® C Economy pack	3 x 150 mL bottles	REF 13065
Applyfix 4 (1:1)		REF 17203

Kettenbach GmbH & Co. KG  
Im Heerfeld 7  
35713 Eschenburg, Germany  
Phone +49 (0)2774 7050  
info@kettenbach.com  
www.kettenbach-dental.com



DEN141114\_3325

**KETTENBACHDENTAL**  
Simply intelligent